

## CSTB antibody

### Product Information

Catalog No.:	FNab02180
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	cystatin B(stefin B)
Synonyms:	Cystatin-B CPI-B Liver thiol proteinase inhibitor Stefin-B CSTB CST6 STFB
Observed MW:	11 kDa
Uniprot ID :	P04080

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab02180(CSTB antibody) at dilution of 1:100



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab02180(CSTB antibody) at dilution of 1:1500