

CYB5R3 antibody

Product Information

Catalog No.:	FNab02186
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	cytochrome b5 type A(microsomal)
Synonyms:	NADH-cytochrome b5 reductase 3 (B5R, Cytochrome b5 reductase) Diaphorase-1 CYB5R3 DIA1
Observed MW:	15 kDa
Uniprot ID :	P00387

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab02186(CYB5A antibody) at dilution of 1:50



human liver tissue were subjected to SDS PAGE followed by western blot with FNab02186(CYB5A antibody) at dilution of 1:800