

CYTH1 antibody

Product Information

Catalog No.:	FNab02190
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cytohesin 1
Synonyms:	Cytohesin-1 PH, SEC7 and coiled-coil domain-containing protein 1 SEC7 homolog B2-1 CYTH1 D17S811E PSCD1
Observed MW:	47 kDa
Uniprot ID :	Q15438

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02190(Cytohesin 1 antibody) at dilution of 1:500