

KRT14 antibody

Product Information


Catalog No.:	FNab02195
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	keratin 14
Synonyms:	Keratin, type I cytoskeletal 14 Cytokeratin-14 (CK-14) Keratin-14 (K14) KRT14
Observed MW:	47-50 kDa
Uniprot ID :	P02533

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP, IHC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	




Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab02195(KRT14 Antibody) at dilution of 1:200



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using FNab02195(KRT14 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



IP Result of anti-Cytokeratin 14 (IP:FNab02195, 3ug; Detection:FNab02195 1:1000) with A431 cells lysate 1200ug.



A431 cells were subjected to SDS PAGE followed by western blot with FNab02195(KRT14 antibody) at dilution of 1:1000