

KRT14 antibody

Product Information

Catalog No.:	FNab02196
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	keratin 14
Synonyms:	Keratin, type I cytoskeletal 14 Cytokeratin-14 (CK-14) Keratin-14 (K14) KRT14
Observed MW:	47-50 kDa
Uniprot ID :	P02533

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab02196(KRT14 Antibody) at dilution of 1:200



A431 cells were subjected to SDS PAGE followed by western blot with FNab02196(KRT14 antibody) at dilution of 1:1000