

## KRT14 antibody

### Product Information


Catalog No.:	FNab02197
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8B5
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	keratin 14
Synonyms:	Keratin, type I cytoskeletal 14 Cytokeratin-14 (CK-14) Keratin-14 (K14) KRT14
Observed MW:	52 kDa
Uniprot ID :	P02533

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:200-1:500; IF: 1:20-1:200
Image:	



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using FNab02197 (Cytokeratin 14 antibody) at dilution of 1:400. heat mediated antigen retrieved with Tris-EDTA buffer (pH9).



A431 cells were subjected to SDS PAGE followed by western blot with FNab02197 (Cytokeratin 14) antibody at dilution of 1:500