

## KRT19-specific antibody

### Product Information

Catalog No.:	FNab02210
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	keratin 19
Synonyms:	Keratin, type I cytoskeletal 19 Cytokeratin-19 (CK-19) Keratin-19 (K19) KRT19
Observed MW:	40 kDa
Uniprot ID :	P08727

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab02210(CK19-specific Antibody) at dilution of 1:200



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02210(KRT19-specific antibody) at dilution of 1:400