

KRT7 antibody

Product Information

Catalog No.:	FNab02217
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	keratin 7
Synonyms:	Keratin, type II cytoskeletal 7 Cytokeratin-7 (CK-7) Keratin-7 (K7) Sarcolectin Type-II keratin Kb7 KRT7 SCL
Observed MW:	54 kDa
Uniprot ID :	P08729

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, FC, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab02217(KRT7 Antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab02217(KRT7 antibody) at dilution of 1:2000