

## KRT7-specific antibody

### Product Information

Catalog No.:	FNab02220
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	keratin 7
Synonyms:	Keratin, type II cytoskeletal 7 Cytokeratin-7 (CK-7) Keratin-7 (K7) Sarcolectin Type-II keratin Kb7 KRT7 SCL
Observed MW:	54 kDa
Uniprot ID :	P08729

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF, IP, IHC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab02220(KRT7-specific antibody) at dilution of 1:100



IP Result of anti-Cytokeratin 7-specific (IP:FNab02220, 4ug; Detection:FNab02220 1:500) with mouse liver tissue lysate 4000ug.



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab02220(KRT7-specific antibody) at dilution of 1:500