

KRT8 antibody

Product Information

Catalog No.:	FNab02222
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	keratin 8
Synonyms:	Keratin, type II cytoskeletal 8 Cytokeratin-8 (CK-8) Keratin-8 (K8) Type-II keratin Kb8 KRT8 CYK8
Observed MW:	54 kDa
Uniprot ID :	P05787

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab02222(KRT8 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab02222(KRT8 Antibody) at dilution of 1:1500