

EIF4G2 antibody

Product Information

Catalog No.:	FNab02239
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	eukaryotic translation initiation factor 4 gamma, 2
Synonyms:	Eukaryotic translation initiation factor 4 gamma 2 (eIF-4-gamma 2, eIF-4G 2, eIF4G 2) Death-associated protein 5 (DAP-5) p97 EIF4G2 DAP5
Observed MW:	97 kDa
Uniprot ID :	P78344

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab02239(DAP5 antibody) at dilution of 1:400