

## DAPK2-Specific antibody

### Product Information

Catalog No.:	FNab02241
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	death-associated protein kinase 2
Synonyms:	Death-associated protein kinase 2 (DAP kinase 2) DAP-kinase-related protein 1 (DRP-1) DAPK2
Observed MW:	43 kDa
Uniprot ID :	Q9UIK4

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:10-1:100
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab02241(DAPK2-Specific antibody) at dilution of 1:200