

DAZ4 antibody

Product Information

Catalog No.:	FNab02250
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


RNA-binding protein that plays an essential role in spermatogenesis. May act by binding to the 3'-UTR of mRNAs and regulating their translation.

Immunogen information

Immunogen:	deleted in azoospermia 4
Synonyms:	Deleted in azoospermia protein 4 DAZ4
Observed MW:	55-60 kDa
Uniprot ID :	Q86SG3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02250(DAZ4 Antibody) at dilution of 1:600