

DAZAP1 antibody

Product Information

Catalog No.:	FNab02251
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


RNA-binding protein, which may be required during spermatogenesis.

Immunogen information

Immunogen:	DAZ associated protein 1
Synonyms:	DAZ-associated protein 1 Deleted in azoospermia-associated protein 1 DAZAP1
Observed MW:	None
Uniprot ID :	Q96EP5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab02251(DAZAP1 antibody) at dilution of 1:50