

BRINP1 antibody

Product Information

Catalog No.:	FNab02253
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	deleted in bladder cancer 1
Synonyms:	BMP/retinoic acid-inducible neural-specific protein 1 Deleted in bladder cancer protein 1 BRINP1 DBC1 DBCCR1 FAM5A
Observed MW:	88 kDa
Uniprot ID :	O60477

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab02253 (DBC1 antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab02253 (DBC1 antibody) at dilution of 1:300