

BRINP1 antibody

Product Information


Catalog No.:	FNab02254
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8C10
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	deleted in bladder cancer 1
Synonyms:	BMP/retinoic acid-inducible neural-specific protein 1 Deleted in bladder cancer protein 1 BRINP1 DBC1 DBCCR1 FAM5A
Observed MW:	88 kDa
Uniprot ID :	O60477

Application

Reactivity:	Human, Pig, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:50-1:500
Image:	



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab02254(DBC1 Antibody) at dilution of 1:1000