

## DCC antibody

### Product Information

Catalog No.:	FNab02262
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	deleted in colorectal carcinoma
Synonyms:	Netrin receptor DCC Colorectal cancer suppressor Immunoglobulin superfamily DCC subclass member 1 Tumor suppressor protein DCC DCC IGDCC1
Observed MW:	180 kDa
Uniprot ID :	P43146

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	

mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab02262( DCC-Specific Antibody) at dilution of 1:600

