

## DCLK1 antibody

### Product Information

Catalog No.:	FNab02266
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probable kinase that may be involved in a calcium-signaling pathway controlling neuronal migration in the developing brain. May also participate in functions of the mature nervous system.

### Immunogen information


Immunogen:	doublecortin-like kinase 1
Synonyms:	Serine/threonine-protein kinase DCLK1 Doublecortin domain-containing protein 3A Doublecortin-like and CAM kinase-like 1 Doublecortin-like kinase 1 DCLK1 DCAMKL1 DCDC3A KIAA0369
Observed MW:	48 kDa, 82 kDa
Uniprot ID :	O15075

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse brain using FNab02266(DCLK1 antibody) at dilution of 1:100



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab02266( DCLK1 Antibody) at dilution of 1:1000