

DCTN6 antibody

Product Information

Catalog No.:	FNab02276
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	dynactin 6
Synonyms:	Dynactin subunit 6 Dynactin subunit p27 Protein WS-3 DCTN6 WS3
Observed MW:	27 kDa
Uniprot ID :	O00399

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab02276(DCTN6 antibody) at dilution of 1:100



IP Result of anti-DCTN6 (IP:FNab02276, 4ug; Detection:FNab02276 1:500) with mouse testis tissue lysate 4800ug.



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab02276(DCTN6 antibody) at dilution of 1:500