

DDAH2 antibody

Product Information


Catalog No.:	FNab02283
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	dimethylarginine dimethylaminohydrolase 2
Synonyms:	Putative hydrolase DDAH2 DDAHII Inactive N(G),N(G)-dimethylarginine dimethylaminohydrolase 2 (DDAH-2, Inactive dimethylarginine dimethylaminohydrolase 2) Protein G6a S-phase protein DDAH2 DDAH G6A NG30
Observed MW:	30 kDa
Uniprot ID :	O95865

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100; IF: 1:10 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab02283(DDAH2 antibody) at dilution of 1:50



Immunofluorescence analysis of HeLa cells using FNab02283(DDAH2 antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab02283(DDAH2 antibody) at dilution of 1:1000