

## DDX25 antibody

### Product Information

Catalog No.:	FNab02302
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


ATP-dependent RNA helicase. Required for mRNA export and translation regulation during spermatid development(By similarity).

### Immunogen information


Immunogen:	DEAD(Asp-Glu-Ala-Asp) box polypeptide 25
Synonyms:	ATP-dependent RNA helicase DDX25 DEAD box protein 25 Gonadotropin-regulated testicular RNA helicase DDX25 GRTH
Observed MW:	50-55 kDa
Uniprot ID :	Q9UHL0

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB : 1:500-1:1000; IP : 1:500-1:1000
Image:	



IP Result of anti-DDX25 (IP: FNab02302, 4ug;  
Detection: FNab02302 1:500) with mouse testis  
tissue lysate 4000ug.



Jurkat cells were subjected to SDS PAGE followed  
by western blot with FNab02302(DDX25  
antibody) at dilution of 1:600