

DDX27 antibody

Product Information

Catalog No.:	FNab02303
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable ATP-dependent RNA helicase. Component of the nucleolar ribosomal RNA(rRNA) processing machinery that regulates 3' end formation of ribosomal 47S rRNA(PubMed:25825154).

Immunogen information

Immunogen:	DEAD(Asp-Glu-Ala-Asp) box polypeptide 27
Synonyms:	Probable ATP-dependent RNA helicase DDX27 DEAD box protein 27 DDX27 cPERP-F RHLP
Observed MW:	90 kDa
Uniprot ID :	Q96GQ7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02303(DDX27 antibody) at dilution of 1:500