

DDX3X antibody

Product Information

Catalog No.:	FNab02305
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	DEAD(Asp-Glu-Ala-Asp) box polypeptide 3, X-linked
Synonyms:	ATP-dependent RNA helicase DDX3X CAP-Rf DEAD box protein 3, X-chromosomal DEAD box, X isoform (DBX) Helicase-like protein 2 (HLP2) DDX3X DBX DDX3
Observed MW:	73 kDa, 60-65 kDa
Uniprot ID :	O00571

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab02305(DDX3 antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab02305(DDX3 Antibody) at dilution of 1:600