

## DEDD2 antibody

### Product Information

Catalog No.:	FNab02324
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a critical role in death receptor-induced apoptosis and may target CASP8 and CASP10 to the nucleus. May regulate degradation of intermediate filaments during apoptosis. May play a role in the general transcription machinery in the nucleus and might be an important regulator of the activity of GTF3C3.

### Immunogen information


Immunogen:	death effector domain containing 2
Synonyms:	DNA-binding death effector domain-containing protein 2 DED-containing protein FLAME-3 FADD-like anti-apoptotic molecule 3 DEDD2 FLAME3
Observed MW:	36 kDa
Uniprot ID :	Q8WXF8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab02324(DED2 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02324(DED2 antibody) at dilution of 1:300