

DEFA1 antibody

Product Information

Catalog No.:	FNab02326
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	defensin, alpha 1
Synonyms:	Neutrophil defensin 1 Defensin, alpha 1 HNP-1 (HP-1, HP1) HP 1-56 Neutrophil defensin 2 Alternative names: HNP-2 (HP-2, HP2) DEFA1 DEF1 DEFA2 MRS
Observed MW:	6 kDa
Uniprot ID :	P59665

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab02326(DEFA1 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02326(DEFA1 Antibody) at dilution of 1:300