

DENND3 antibody

Product Information

Catalog No.:	FNab02333
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Guanine nucleotide exchange factor(GEF) activating RAB12. Promotes the exchange of GDP to GTP, converting inactive GDP-bound RAB12 into its active GTP-bound form. Thereby, may play a role in protein transport from recycling endosomes to lysosomes regulating, for instance, the degradation of the transferrin receptor.

Immunogen information

Immunogen:	DENN/MADD domain containing 3
Synonyms:	DENN domain-containing protein 3 DENND3 KIAA0870
Observed MW:	136 kDa
Uniprot ID :	A2RUS2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	



mouse stomach tissue were subjected to SDS
PAGE followed by western blot with
FNab02333(DENND3 Antibody) at dilution of
1:1000