
DES antibody

Product Information


Catalog No.:	FNab02342
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6H3
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	desmin
Synonyms:	Desmin DES
Observed MW:	50-55 kDa
Uniprot ID :	P17661

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab02342(DES antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab02342(DES antibody) at dilution of 1:1000