

DSC2 antibody

Product Information

Catalog No.:	FNab02344
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7D7
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	desmocollin 2
Synonyms:	Desmocollin-2 Cadherin family member 2 Desmocollin-3 Desmosomal glycoprotein II Desmosomal glycoprotein III DSC2 CDHF2 DSC3
Observed MW:	110 kDa
Uniprot ID :	Q02487

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:1000-1:4000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02344(Desmocollin 2 antibody) at dilution of 1:2000