

## DGCR2 antibody

### Product Information

Catalog No.:	FNab02350
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Putative adhesion receptor, that could be involved in cell-cell or cell-matrix interactions required for normal cell differentiation and migration.

### Immunogen information

Immunogen:	DiGeorge syndrome critical region gene 2
Synonyms:	Integral membrane protein DGCR2/IDD DGCR2 IDD KIAA0163
Observed MW:	61-65 kDa
Uniprot ID :	P98153

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:200