

DGUOK antibody

Product Information

Catalog No.:	FNab02354
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	deoxyguanosine kinase
Synonyms:	Deoxyguanosine kinase, mitochondrial Deoxyadenosine kinase, mitochondrial DGUOK DGK
Observed MW:	28 kDa
Uniprot ID :	Q16854

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab02354(DGUOK antibody) at dilution of 1:100



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab02354(DGUOK antibody) at dilution of 1:300