

DGUOK antibody

Product Information

Catalog No.:	FNab02355
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	deoxyguanosine kinase
Synonyms:	Deoxyguanosine kinase, mitochondrial Deoxyadenosine kinase, mitochondrial DGUOK DGK
Observed MW:	32 kDa
Uniprot ID :	Q16854

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IP
Recommended dilution:IP:	1:200-1:2000
Image:	



IP Result of anti-dGK (IP:FNab02355, 3ug;
Detection:FNab02355 1:300) with mouse kidney
tissue lysate 4000ug.