

DGKB antibody

Product Information

Catalog No.:	FNab02357
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Exhibits high phosphorylation activity for long-chain diacylglycerols.

Immunogen information

Immunogen:	diacylglycerol kinase, beta 90kDa
Synonyms:	Diacylglycerol kinase beta (DAG kinase beta) 90 kDa diacylglycerol kinase Diglyceride kinase beta (DGK-beta) DGKB DAGK2 KIAA0718
Observed MW:	91 kDa
Uniprot ID :	Q9Y6T7

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02357(DGKB Antibody) at dilution of 1:1000