

## DHDH antibody

### Product Information

Catalog No.:	FNab02364
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	dihydrodiol dehydrogenase(dimeric)
Synonyms:	Trans-1,2-dihydrobenzene-1,2-diol dehydrogenase D-xylose 1-dehydrogenase D-xylose-NADP dehydrogenase Dimeric dihydrodiol dehydrogenase Hum2DD DHDH 2DD
Observed MW:	39 kDa
Uniprot ID :	Q9UQ10

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab02364(DHDH antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02364(DHDH antibody) at dilution of 1:500