

## DHH antibody

### Product Information

Catalog No.:	FNab02366
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Intercellular signal essential for a variety of patterning events during development. May function as a spermatocyte survival factor in the testes. Essential for testes development.

### Immunogen information

Immunogen:	desert hedgehog homolog(Drosophila)
Synonyms:	Desert hedgehog protein (DHH) HHG-3 Desert hedgehog protein N-product (DHH-N) DHH
Observed MW:	42 kDa
Uniprot ID :	O43323

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab02366(DHH antibody) at dilution of 1:1000