

## DHX15 antibody

### Product Information

Catalog No.:	FNab02375
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Pre-mRNA processing factor involved in disassembly of spliceosomes after the release of mature mRNA. In cooperation with TFIP11 seem to be involved in the transition of the U2, U5 and U6 snRNP-containing IL complex to the snRNP-free IS complex leading to efficient debranching and turnover of excised introns.

### Immunogen information

Immunogen:	DEAH(Asp-Glu-Ala-His) box polypeptide 15
Synonyms:	ATP-dependent RNA helicase DHX15 ATP-dependent RNA helicase #46 DEAH box protein 15 Splicing factor Prp43 (hPrp43) DHX15 DBP1 DDX15
Observed MW:	91 kDa
Uniprot ID :	O43143

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab02375( DHX15 Antibody) at dilution of 1:600