

DHX32 antibody

Product Information

Catalog No.:	FNab02378
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	DEAH(Asp-Glu-Ala-His) box polypeptide 32
Synonyms:	Putative pre-mRNA-splicing factor ATP-dependent RNA helicase DHX32 DEAD/H box 32 DEAD/H helicase-like protein 1 (DHLPI) DEAH box protein 32 HuDDX32 DHX32 DDX32
Observed MW:	84 kDa, 75 kDa
Uniprot ID :	Q7L7V1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab02378(DHX32 antibody) at dilution of 1:100



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab02378(DHX32 antibody) at dilution of 1:800