

DIABLO antibody

Product Information


Catalog No.:	FNab02384
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	diablo homolog(Drosophila)
Synonyms:	Diablo IAP-binding mitochondrial protein Diablo homolog, mitochondrial Direct IAP-binding protein with low pI Second mitochondria-derived activator of caspase (Smac) Diablo IAP-binding mitochondrial protein, cleaved form DIABLO SMAC
Observed MW:	20 kDa
Uniprot ID :	Q9NR28

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab02384(DIABLO antibody) at dilution of 1:50



IP Result of anti-DIABLO (IP:FNab02384, 4ug; Detection:FNab02384 1:500) with HEK-293 cells lysate 2000ug.



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02384(DIABLO antibody) at dilution of 1:1000