

DIDO1 antibody

Product Information

Catalog No.:	FNab02389
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Putative transcription factor, weakly pro-apoptotic when overexpressed(By similarity). Tumor suppressor. Required for early embryonic stem cell development. Isoform 2: Displaces isoform 4 at the onset of differentiation, required for repression of stemness genes.

Immunogen information

Immunogen:	death inducer-obliterator 1
Synonyms:	Death-inducer obliterator 1 (DIO-1, hDido1) Death-associated transcription factor 1 (DATF-1) DIDO1 C20orf158 DATF1 KIAA0333
Observed MW:	None
Uniprot ID :	Q9BTC0

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab02389(DIDO1 antibody) at dilution of 1:3000