

## DIXDC1 antibody

### Product Information

Catalog No.:	FNab02399
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Positive effector of the Wnt signaling pathway; activates WNT3A signaling via DVL2. Regulates JNK activation by AXIN1 and DVL2.

### Immunogen information


Immunogen:	DIX domain containing 1
Synonyms:	Dixin Coiled-coil protein DIX1 (Coiled-coil-DIX1) DIX domain-containing protein 1 DIXDC1 CCD1 KIAA1735
Observed MW:	70-75 kDa
Uniprot ID :	Q155Q3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab02399(DIXDC1 antibody) at dilution of 1:100



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab02399(DIXDC1 antibody) at dilution of 1:1500