

DLG5 antibody

Product Information

Catalog No.:	FNab02411
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role at the plasma membrane in the maintenance of the structure of epithelial cells and in the transmission of extracellular signals to the membrane and cytoskeleton.

Immunogen information

Immunogen:	discs, large homolog 5(Drosophila)
Synonyms:	Disks large homolog 5 Discs large protein P-dlg Placenta and prostate DLG DLG5 KIAA0583 PDLG
Observed MW:	75 kDa, 105 kDa
Uniprot ID :	Q8TDM6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab02411(DLG5 antibody) at dilution of 1:500