

## DLGAP3 antibody

### Product Information

Catalog No.:	FNab02412
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in the molecular organization of synapses and neuronal cell signaling. Could be an adapter protein linking ion channel to the subsynaptic cytoskeleton. May induce enrichment of PSD-95/SAP90 at the plasma membrane.

### Immunogen information

Immunogen:	discs, large(Drosophila) homolog-associated protein 3
Synonyms:	Disks large-associated protein 3 (DAP-3) PSD-95/SAP90-binding protein 3 SAP90/PSD-95-associated protein 3 (SAPAP3) DLGAP3 DAP3
Observed MW:	106-120 kDa
Uniprot ID :	O95886

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02412(DLGAP3 antibody) at dilution of 1:500