

## DLK2 antibody

### Product Information

Catalog No.:	FNab02414
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Regulates adipogenesis.

### Immunogen information

Immunogen:	delta-like 2 homolog(Drosophila)
Synonyms:	Protein delta homolog 2 (DLK-2) Epidermal growth factor-like protein 9 (EGF-like protein 9) DLK2 EGFL9
Observed MW:	41 kDa
Uniprot ID :	Q6UY11

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab02414(DLK2 antibody) at dilution of 1:500