

DLL1-Specific antibody

Product Information

Catalog No.:	FNab02415
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	delta-like 1(Drosophila)
Synonyms:	Delta-like protein 1 Drosophila Delta homolog 1 (Delta1, H-Delta-1) DLL1
Observed MW:	80 kDa
Uniprot ID :	O00548

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab02415(DLL1-Specific antibody) at dilution of 1:500