

DLX3 antibody

Product Information

Catalog No.:	FNab02418
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.

Immunogen information


Immunogen:	distal-less homeobox 3
Synonyms:	Homeobox protein DLX-3 DLX3
Observed MW:	32 kDa
Uniprot ID :	O60479

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNa02418(DLX3 antibody) at dilution of 1:100



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNa02418(DLX3 antibody) at dilution of 1:1000