

DNAH6 antibody

Product Information

Catalog No.:	FNab02436
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Force generating protein of respiratory cilia. Produces force towards the minus ends of microtubules. Dynein has ATPase activity; the force-producing power stroke is thought to occur on release of ADP(By similarity).

Immunogen information

Immunogen:	dynein, axonemal, heavy chain 6
Synonyms:	Dynein axonemal heavy chain 6 Axonemal beta dynein heavy chain 6 Ciliary dynein heavy chain 6 DNAH6 DNAHC6 DNHL1 HL2 KIAA1697
Observed MW:	None
Uniprot ID :	Q9C0G6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab02436(DNAH6 Antibody) at dilution of 1:50