

DNAJB5 antibody

Product Information

Catalog No.:	FNab02451
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	DnaJ(Hsp40) homolog, subfamily B, member 5
Synonyms:	DnaJ homolog subfamily B member 5 Heat shock protein Hsp40-2 Heat shock protein Hsp40-3 Heat shock protein cognate 40 (Hsc40) DNAJB5 HSC40
Observed MW:	39 kDa
Uniprot ID :	O75953

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab02451(DNAJB5 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02451(DNAJB5 antibody) at dilution of 1:500