

DNAJC11 antibody

Product Information

Catalog No.:	FNab02458
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	DnaJ(Hsp40) homolog, subfamily C, member 11
Synonyms:	DnaJ homolog subfamily C member 11 DNAJC11
Observed MW:	63 kDa
Uniprot ID :	Q9NVH1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using FNab02458(DNAJC11 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02458(DNAJC11 antibody) at dilution of 1:1000