

DNAJC12 antibody

Product Information

Catalog No.:	FNab02459
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	DnaJ(Hsp40) homolog, subfamily C, member 12
Synonyms:	DnaJ homolog subfamily C member 12 J domain-containing protein 1 DNAJC12 JDP1
Observed MW:	23 kDa
Uniprot ID :	Q9UKB3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab02459(DNAJC12 antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab02459(DNAJC12 antibody) at dilution of 1:300