

DNAJC2 antibody

Product Information

Catalog No.:	FNab02463
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	DnaJ(Hsp40) homolog, subfamily C, member 2
Synonyms:	DnaJ homolog subfamily C member 2 M-phase phosphoprotein 11 Zuotin-related factor 1 DnaJ homolog subfamily C member 2, N-terminally processed DNAJC2 MPHOSPH11 MPP11 ZRF1
Observed MW:	66-75 kDa
Uniprot ID :	Q99543

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab02463(DNAJC2/MPP11 antibody) at dilution of 1:300