

DNAJC21 antibody

Product Information

Catalog No.:	FNab02464
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as a co-chaperone for HSP70.

Immunogen information


Immunogen:	DnaJ(Hsp40) homolog, subfamily C, member 21
Synonyms:	DnaJ homolog subfamily C member 21 DnaJ homolog subfamily A member 5 Protein GS3 DNAJC21 DNAJA5
Observed MW:	70 kDa
Uniprot ID :	Q5F1R6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab02464(DNAJC21 Antibody) at dilution of 1:50



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab02464(DNAJC21 antibody) at dilution of 1:2000